

## **Preparation of B-Lymphocyte Lysates for Cyclic AMP Determination**

AfCS Procedure Protocol ID PP0000001200

Version 1, 03/07/02

The method chosen for measuring the content of cyclic adenosine 3',5'-monophosphate (cyclic AMP or cAMP) in splenic B lymphocytes (B cells) is an enzyme-linked immunoassay system developed by Amersham Biosciences. An overview of the assay, [cAMP Biotrak EIA](#), and its instructions, [Biotrak Protocol](#), are available through the indicated links. The procedure described below provides a sufficient sample for several determinations of cAMP using the acetylation protocol.

### **Treatment of Cells and Preparation of Extracts**

1. Suspend freshly isolated splenic B cells at  $16.7 \times 10^6$  cells/ml in Supplemented Iscove's Modified Dulbecco's Medium (SIMDM) and distribute 0.090 ml into individual wells of 96-well tissue culture plates as needed.
2. Incubate at 37 °C in air and 5% CO<sub>2</sub> for 1 hr.
3. Transfer the plate to an environmental chamber containing air at 37 °C.
4. At timed intervals, add 0.010 ml of ligand working stock (10X final concentration in SIMDM) or vehicle control working stock (in SIMDM) to appropriate wells to begin treatments. (Note: vehicle controls constitute matching dilutions of solvents in which ligands are dissolved and stored.)
5. Mix with a microtiter plate shaker for 2 sec at 350 rpm.
6. Incubate at 37 °C in the environmental chamber.
7. At the desired times, end treatments with the addition of 186 µl of 100% ethanol and mix by pipetting up and down 3 times. (A multichannel pipette is used for parallel time courses with several ligands.)
8. After all of the samples on the plate have been stopped, allow lysates to sit on ice for 5 min. Centrifuge the plate at 2000 x g for 15 min at 4°C. Transfer the lysates to barcoded, 0.6-ml Eppendorf tubes.
9. Add 100 µl of 65% ethanol to the same wells and mix by tilting the plate in a circular motion to rinse the wells. Centrifuge the plate at 2000 x g for 5 min at 4°C.
10. Combine the supernatant (step 9) with the appropriate lysates (step 8) in the 0.6 ml Eppendorf tubes.
11. Dry the samples in a vacuum with centrifugation and heat (Speed Vac) until the samples are completely dry (about 2 to 3 hr).
12. Secure caps on the samples and store at -80 °C until determination of cAMP by enzyme immunoassay (EIA).

### **Reagents and Materials**

Supplemented Iscove's Modified Dulbecco's Medium (SIMDM): AfCS Solution Protocol ID PS0000005600

96-well, tissue culture plates: Falcon; catalog no. 353072

Ethanol, 100%: Aaper Alcohol and Chemical Co.; catalog no. 030801

Eppendorf tubes, 0.6 ml: Sarstedt; catalog no. 72.699

Ethanol, 65%: AfCS Solution Protocol ID PS0000007500

Speed Vac System: Savant; catalog no. SS21

**Author:** Keng-Mean Lin

**Date:** 03/8/02

**Approved:** Paul C. Sternweis