

## Quantification of Apoptosis and the Cell Cycle Distribution of Primary B Cells Using Propidium Iodide

AfCS Procedure Protocol ID PP0000002400  
Version 1, 02/13/02

The water soluble, DNA intercalator, propidium iodide (PI), is used to bind to DNA after permeabilization of cells with NP40. The amount of dye bound correlates with the content of DNA within a given cell. Once cells are stained, they are analyzed on a flow cytometer. The relative content of DNA indicates the distribution of a population of cells throughout the cell cycle. For example, cells in the  $G_0/G_1$  phases of the cell cycle are diploid, or have a DNA content of  $2n$ . Cells within the  $G_2/M$  phases have a DNA content of  $4n$ , while S-phase cells have a DNA content greater than  $2n$  and less than  $4n$ . Cells that are sub-diploid ( $<2n$ ) are apoptotic and can also be quantified using this method.

### Procedures

1. Remove  $1 \times 10^6$  cells from the appropriate samples and place into microfuge tubes.
2. Centrifuge in a microfuge at  $400 \times g$  for 5 min at  $4^\circ\text{C}$ .
3. Aspirate the supernatants.
4. Cap the tubes. Invert the tubes (approximately 150 degrees) and flick the conical end with a finger 3 to 4 times to loosen the pellets.
5. Add  $500 \mu\text{l}$  of fluorescence-activated cell sorting buffer (FACS buffer) to each sample and gently invert approximately 3 times to wash the cells.
6. Centrifuge in a microfuge for 5 min at  $400 \times g$  at  $4^\circ\text{C}$ .
7. Aspirate the supernatants.
8. Add  $250 \mu\text{l}$  of propidium iodide staining solution (PI solution) to each sample.
9. Place the tubes on ice in an ice bucket and cover.
10. Incubate for 30 min.
11. Add  $500 \mu\text{l}$  of FACS buffer to each tube.
12. Transfer each sample from a microfuge tube to a FACS tube.
13. Mix and analyze by passage through a FACSCalibur flow cytometer.

### Reagents and Materials

Fluorescence-activated cell sorting buffer (FACS buffer): AfCS Solution Protocol ID PS0000002900

Propidium iodide staining solution (PI solution): AfCS Solution Protocol ID PP0000004400

FACS tubes (Falcon polystyrene round-bottom tubes, 12 x 75 mm): Falcon; catalog no. 35-2054

FACSCalibur flow cytometer: Becton Dickinson; catalog no. 343023

**Author:** Robert Hsueh/Tamara Roach

**Date:** 02/15/02

**Approved:** Paul Sternweis/Bill Seaman